

INSTALLATION INSTRUCTIONS *for Additional Kit Purchases*

Follow these instructions to install newly purchased kits to a machine where any version of VHI PC-Kits has been previously installed.

After installation, you may be required to activate PC-Kits before it can be used. Activation is available via the Internet 24 hours a day or by phone during the support hours listed below.

1. Exit all running software and temporarily disable any virus protection software you may be using.
2. Place the VHI PC-Kits CD into your CD drive. The setup program will automatically start. If the setup program does not start, go to **My Computer** and right-click on the CD drive. Choose **AutoPlay** from the menu.
3. From the splash screen that appears, click on **Install or Remove PC-Kits** (see Figure 1).
4. Depending on which version of VHI PC-Kits is currently installed on your machine, an upgrade program may now run. Let the upgrade complete and restart your computer if requested to do so. After restarting, repeat steps 2 and 3, then proceed to step 5.
5. Follow the on screen instructions until asked whether you want to **Remove, Repair** or **Modify** (see Figure 2). Choose **Modify** and click **Next**.
6. Continue to follow the on screen instructions until the **Serial Numbers** dialog appears (see Figure 3). Enter the new serial numbers that you received with your new kit purchase. After entering each serial number, press the **Add** button.
7. After you have added your new serial number(s), click **Next** and continue to follow the remaining on screen instructions. If asked to reboot at the end of setup, then do so before continuing to step 8.

Directions continue on opposite side.




Figure 1. Splash Screen




Figure 2. Remove, Repair or Modify Dialog




Figure 3. The Serial Numbers Dialog

8. Start VHI PC-Kits by double clicking the desktop icon. If you see the dialog shown in Figure 4, you will need to activate your PC-Kits installation before you can use it.
9. If you are connected to the Internet, select ***Yes, I want to Activate over the Internet.*** If you are not connected to the Internet, select ***Yes, I want to Activate by telephone.*** Click **OK**.
10. Follow the on screen instructions to perform the activation process. You only need to activate PC-Kits once each time you purchase new kits.
11. If you received a new CD with this purchase please discard all your old VHI PC-Kits CDs. It is only necessary to keep the latest CD. Be sure to keep this instruction sheet and all serial numbers in a safe place with your PC-Kits box or manual.




Figure 4. The Activation Dialog

**If you encounter any difficulties
during the installation process,
please call Health Promotion Resources
on 03 9888 1829
for product support.
Monday - Friday 9.00am - 4.00pm AEST**